

REF			SYSTEM
08946728190	08946728500	100	cobas e 402 cobas e 801

## English

### System information

Short name	ACN (application code number)
ACTH	10206

### Intended use

Immunoassay for the in vitro quantitative determination of adrenocorticotrophic hormone (ACTH) in human EDTA plasma.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

### Summary

Adrenocorticotrophic hormone or corticotropin is a peptide hormone consisting of 39 amino acids. It is produced in the anterior pituitary of the brain as part of the precursor molecule pro-opiomelanocortin (POMC). Tissue-specific cleavage results in ACTH and a range of related peptides.<sup>1,2</sup> ACTH stimulates formation and secretion of glucocorticoids (especially cortisol) by the adrenal cortex.

The glucocorticoid production is regulated by various factors.<sup>3,4,5,6</sup> After stimulation (e.g. by physical effort or by the internal body clock), the hypothalamus secretes CRH (corticotropin releasing hormone). CRH acts on the pituitary, which in turn synthesizes and secretes ACTH. Finally, ACTH stimulates secretion of the glucocorticoids by the adrenals. High concentrations of glucocorticoids in the blood inhibit secretion of CRH and ACTH via a negative feedback mechanism.

ACTH concentrations show a diurnal variation with high levels in the morning and low levels in the evening. Therefore, as with cortisol, it is important to know the collection time of the plasma sample for interpretation of the results.

Plasma ACTH measurements are useful in the differential diagnosis of Cushing's disease (ACTH hypersecretion), autonomous ACTH producing pituitary tissue (e.g. Nelson's syndrome), hypopituitarism with ACTH deficiency and ectopic ACTH syndrome.<sup>7,8</sup> In addition to cortisol measurements, ACTH determinations can be used together with suppression or stimulation tests to diagnose the origin of glucocorticoid overproduction. Similarly, ACTH measurements can be employed to facilitate differential diagnosis of adrenocortical insufficiency (Addison's disease).<sup>9</sup>

ACTH not produced by the pituitary gland is known as ectopic ACTH;<sup>10</sup> this is often associated with small cell carcinoma of the lung. In rare cases ectopic ACTH can be caused by thymic tumors, pancreatic adenocarcinomas, or bronchial carcinoids. These tumors often secrete ACTH precursors (POMC and pro-ACTH).

The Elecsys ACTH assay employs 2 monoclonal antibodies specific for ACTH (9-12) and for the C-terminal region (ACTH 36-39).

Due to common antigenic structure, the antibodies recognize intact biologically active ACTH 1-39 and the ACTH precursors POMC and pro-ACTH.<sup>2</sup>

### Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: 30 µL of sample, a biotinylated monoclonal ACTH-specific antibody, and a monoclonal ACTH-specific antibody labeled with a ruthenium complex<sup>9</sup> react to form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the **cobas** link.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)<sub>3</sub><sup>2+</sup>)

### Reagents - working solutions

The **cobas e** pack is labeled as ACTH.

- M Streptavidin-coated microparticles, 1 bottle, 5.8 mL:  
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-ACTH-Ab~biotin, 1 bottle, 7.2 mL:  
Biotinylated monoclonal anti-ACTH antibody (mouse) 0.3 mg/L;  
MES<sup>b)</sup> buffer 50 mmol/L, pH 6.2; preservative.
- R2 Anti-ACTH-Ab~Ru(bpy)<sub>3</sub><sup>2+</sup>, 1 bottle, 7.2 mL:  
Monoclonal anti-ACTH antibody (mouse) labeled with ruthenium complex 0.3 mg/L; MES buffer 50 mmol/L, pH 6.2; preservative.

b) MES = 2-morpholino-ethane sulfonic acid

### Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



### Warning

H317 May cause an allergic skin reaction.

### Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

### Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

### Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

### Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is available via the **cobas** link.

### Storage and stability

Store at 2-8 °C.





Cross-reactant	Concentration of cross-reactant pg/mL	Apparent ACTH pg/mL	Change in ACTH concentration pg/mL	Cross-reactivity %
None; reference	0	55.4	not applicable	not applicable
ACTH 22-39	50000	7.58	-47.8	-0.096
	5000	37.5	-17.9	-0.357
	500	52.9	-2.5	-0.491
ACTH 1-13 (alpha-MSH)	50000	29.2	-26.2	-0.052
	5000	51.4	-4.0	-0.080
	500	55.3	-0.1	-0.022

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2021, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

+800 5505 6606



## References

- 1 Cone RD. Anatomy and regulation of the central melanocortin system. *Nature Neurosci* 2005;8:571-578.
- 2 Talbot JA, Kane JW, White A. Analytical and clinical aspects of adrenocorticotrophin determination. *Ann Clin Biochem* 2003;40:453-471.
- 3 Jacobson L. Hypothalamic-pituitary-adrenocortical axis regulation. *Endocrinol Metab Clin North Am* 2005;34:271-292.
- 4 Arit W, Stewart PM. Adrenal corticosteroid biosynthesis, metabolism, and action. *Endocrinol Metab Clin North Am* 2005;34:293-313.
- 5 Lin L, Achermann JC. The Adrenal. *Horm Res* 2004;62:22-29.
- 6 Engelmann M, Landgraf R, Wotjak CT. The hypothalamic-neurohypophysial system regulates the hypothalamic-pituitary-adrenal axis under stress: an old concept revisited. *Front Neuroendocrinol* 2004;25:132-149.
- 7 Beauregard C, Dickstein G, Lacroix A. Classic and recent etiologies of Cushing's syndrome: diagnosis and therapy. *Treat Endocrinol* 2002;1:79-94.
- 8 Lindsay JR, Nieman LK. Differential diagnosis and imaging in Cushing's syndrome. *Endocrinol Metab Clin North Am* 2005;34:403-421.
- 9 Napier C, Pearce SHS. Current and emerging therapies for Addison's disease. *Curr Opin Endocrinol Diabetes Obes* 2014;21:147-153.
- 10 Oliver RL, Davis JR, White A. Characterization of ACTH related peptides in ectopic Cushing's syndrome. *Pituitary* 2003;6:119-126.
- 11 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [dialog.roche.com](http://dialog.roche.com) for definition of symbols used):

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number